

RiboFlow[®]

BIOMOLECULAR
PATHOGEN DETECTION

**AWARD WINNING
TECHNOLOGY**



- Reliable
- Simple
- Rapid
- Safe
- Economic

MOLECULAR FLOW BASED PATHOGEN DETECTION

SY-LAB's RiboFlow® technology signals relevant bacterial contamination on-site – early, quickly, and reliably reducing the need for sophisticated sampling and amplification techniques.

Rapid results are critical when testing for harmful pathogens in foods, toiletries and cosmetics. RiboFlow® technology represents a new, sensitive and user-friendly class of diagnostics:

Molecular Flow-Based Detection. The unique horizontal flow technology allows for the specificity of PCR tests, but with the ease of use associated with lateral flow assays. The result is an effective and rapid system that minimizes sample preparation, speeds time to results, and

provides easy to interpret data for the end user.

The RiboFlow® technology is an advance over current diagnostic capabilities, providing molecular results in an easy to use format and without sophisticated sample preparation steps.

NO gels, NO expensive instrumentation, NO inconclusive and difficult to interpret results.

Fast, easier, and more affordable than conventional molecular tests – SY-LAB's RiboFlow® Technology is the next generation in powerful diagnostic tools helping manufacturers and government regulators to keep the global food and toiletry supply safe.



simple testing



reduced workload



rapid process intervention

YOUR BENEFITS

- **RELIABLE**

- proven biomolecular technology
- novel rRNA based highly specific pathogen detection system

- **SIMPLE**

- simple to perform
- no sophisticated instrumentation necessary

- **RAPID**

- results available in < 30 minutes

- **SAFE**

- unsurpassed specificity
- only viable microbes detected

- **ECONOMIC**

- no capital costs
- minimal hands on time
- affordable assay costs
- stable at ambient temperature



simple documentation



rapid results



earlier product release

TEST PROCEDURE

1) CENTRIFUGATION

(only for liquid enrichments)

2) LYSIS, DENATURATION AND NEUTRALISATION

of released nucleic acids

3) SAMPLE APPLICATION AND INCUBATION

transfer lysate to a RiboFlow test device
and incubate at 48°C for 15 minutes

4) RESULT

visual evaluation and photographic documentation



Fig:
negative test result (l.)
positive test result (r.)

AVAILABLE

RiboFlow® TEST SYSTEMS

RiboFlow® – Salmonella

RiboFlow® – Cronobacter

RiboFlow® – Listeria Twin
(Listeria spp. and
L. monocytogenes
combitest)

RiboFlow® – Campylobacter

RiboFlow® – Vibrio parahaemolyticus

RiboFlow® – Staphylococcus aureus

RiboFlow® – Enterobacteriaceae

RiboFlow® – E.coli

RiboFlow® – Pseudomonas aeruginosa

RiboFlow® – Candida albicans

RiboFlow® – Cosmetics

RiboFlow® – P/En (Pseudomonas spp.
and Enterobacteriaceae
combitest)

RiboFlow® – Ali Twin (Alicyclobacillus
acidoterrestris and Alicyclo-
bacillus spp. combitest)



STARTING MATERIAL

- **Enrichment cultures according to product specification**
- **Isolated colonies from (selective-) agar plates**

RiboFlow®

ACCESSORIES

- **Mini-Incubator IL 10**
- **Thermo-Manipulation Plate**
- **Microcentrifuge**



SY-LAB

WE DEVELOP **SOLUTIONS**
FOR A **SAFER FUTURE**

SY-LAB Geräte GmbH
Tullnerbachstraße 61-65
3011 Neupurkersdorf, Austria
www.sylab.com . sales@sylab.com

Version 1.1 gb